


ProTect™ COVID-19 PCR Kit

Accurate detection of the Sars-CoV-2 virus

Managing COVID-19: testing is still the key

The coronavirus disease 2019 (COVID-19) is an infectious disease of the respiratory system that is caused by the SARS-CoV-2 virus, which is closely related to the SARS virus. The outbreak started in Wuhan at the end of 2019 and to date (16 March 2020), has infected about 170,000 people worldwide and resulted in more than 6,000 deaths. The epicenter of the outbreak has also shifted from Hubei, China to Europe. WHO has also officially declared the outbreak a pandemic. To contain the outbreak, it is essential to diagnose infected patients so that they can be isolated and treated. The ability to gain access to affordable and good quality diagnostic kits is critical.



The image is a screenshot of a BBC News website page. At the top, there is a navigation bar with the BBC logo, a 'Sign in' button, and various category links: News, Sport, Reel, Worklife, Travel, Future, Culture, Music, TV, Weather, and Sounds. Below this is a red banner with the word 'NEWS' in white. Underneath the banner is a secondary navigation bar with links to Home, Video, World, Asia, UK, Business, Tech, Science, Stories, Entertainment & Arts, Health, World News TV, In Pictures, Reality Check, and Newsbeat. Below this is a third navigation bar with links to World, Africa, Australia, Europe, Latin America, Middle East, and US & Canada. The main content area features a headline: 'WHO head: 'Our key message is: test, test, test''. Below the headline is a sub-headline: 'World Health Organisation head Tedros Adhanom Ghebreyesus says there has not been an urgent enough escalation in testing, isolation and contact tracing, which should be the "backbone" of the global response.' Below the sub-headline is a paragraph: 'He said it is not possible to "fight a fire blindfolded", and social distancing measures and handwashing will not alone extinguish the epidemic.' To the right of the text is a video player showing a man in a suit and glasses speaking. The video player has a title bar that says 'test, test, test. Test' and a subtitle that says 'Dr Tedros Adhanom Ghebreyesus WHO DIRECTOR-GENERAL'. The date '16 Mar 2020' is displayed at the bottom left of the article.

WHO head: 'Our key message is: test, test, test'

World Health Organisation head Tedros Adhanom Ghebreyesus says there has not been an urgent enough escalation in testing, isolation and contact tracing, which should be the "backbone" of the global response.

He said it is not possible to "fight a fire blindfolded", and social distancing measures and handwashing will not alone extinguish the epidemic.


16 Mar 2020

PCR is still the gold standard for COVID-19 testing

What is PCR?

It stands for polymerase chain reaction. It is a method to identify the presence of the virus. It “searches” for a particular gene unique to the virus. Once that gene is found, it then makes many copies of that gene (much like a photocopying machine). Eventually, the process makes millions of copies of that gene, which then generates sufficient fluorescent signals to be detected.

PCR can detect the presence of the virus once the person is infected. Hence, it is still the most accurate method for COVID-19 infection.



About JN Medsys



JN Medsys is a Singapore based manufacturer of high-quality PCR and digital PCR products which comprises instruments, consumables and assay kits. It has ISO13485 Quality Management System that allows it to research, develop, manufacture and commercialise the products for in-vitro diagnostics purposes.

THE STRAITS TIMES

B16 | MONEY

WEDNESDAY, OCTOBER 15, 2020

Faster DNA analysis with S'pore product

Local firm to launch its system in the US to compete with industry bigwigs

By RACHAEL BOON

A HOME-GROWN biomedical instrumentation company is launching its first product, which it claims can analyse DNA much more quickly than rivals.
The product, which will be used in laboratories for research, is unveiled on a global stage. The firm, JN Medsys, opened a system that performs known as digital polymerase chain reaction (PCR) to detect deoxyribonucleic acid.

doing the same thing. All are US companies, including current market leader, Bio-Rad laboratories.
JN Medsys' system, which costs about US\$40,000 (S\$51,000), is faster as the entire process takes about 2 1/2 hours compared with up to six hours in others, he said.
Dr Ng added: "JN Medsys' system is easier to use and can analyse up to three times more samples at half the price."

Digital PCR can be used to



JN Medsys chief executive Johanson Ng hopes to start shipping the Clarity digital PCR system by next year. PHOTO: DIOS VINCOY JR FOR THE STRAITS TIMES

working with Dr Tan Min-H principal investigator at A*Star Institute of Bioengineering & Nanotechnology and visiting

ASIANSCIENTIST

TOP NEWS IN THE LAB HEALTH TECHNOLOGY PHARMA ACADEMIA FEATURES

Developing Diagnostics In Double-Quick Time

JN Medsys has developed two new coronavirus diagnostic kits to meet the surge in demand for COVID-19 diagnostics.

SHARE



AsianScientist (Mar. 20, 2020) – All over the world, countries are sinking into a public health crisis as the COVID-19 outbreak develops into a global pandemic. Highly infectious, the coronavirus SARS-CoV-2 can spread from person to person even before symptoms appear. Accurate, widespread diagnosis is therefore critical to curbing new viral outbreaks, allowing infected individuals to be identified and quickly quarantined to stop the train of transmission.

South Korea, one of the first countries to experience the epidemic outside China, has shown that widespread testing might be a good solution to contain the virus. The country's elaborate testing program, which saw more than 270,000 people tested thus far, has helped to slow down the influx of new cases since its record high of 909 cases on February 29, 2020.

JN Medsys Special

DEMOCRATIZING THE DIGITAL REVOLUTION



What is dPCR and how does it work?

PCR is what is known as a PCR amplification method. That's where digital PCR comes in.
Another group of people who would benefit are those working on infectious diseases. Say for example, you are studying pathogen like hepatitis B or hepatitis C, where you only need to be able to quantify the viral load. dPCR will be very helpful. There are just lots of other applications.

How does JN Medsys differ from other products on the market?

JN Medsys' digital PCR products are designed to be simple, easy to use and easy to integrate into existing workflows. For example, one of the key reasons we decided to go digital was to make it easy to use in a lab setting.

What are the advantages of

How much lower and more affordable are JN Medsys products?

JN Medsys' digital PCR products are much more affordable than other products on the market. In fact, our products are up to 10 times cheaper than other products on the market. Our products are also much more reliable and accurate than other products on the market.

What are your near-term plans for the company?

JN Medsys' digital PCR products are designed to be simple, easy to use and easy to integrate into existing workflows. For example, one of the key reasons we decided to go digital was to make it easy to use in a lab setting.

JN Medsys' digital PCR products are designed to be simple, easy to use and easy to integrate into existing workflows. For example, one of the key reasons we decided to go digital was to make it easy to use in a lab setting.



"Our overwhelming focus right now is to build a good product."
Dr. Johanson Ng
Chief Executive Officer
JN Medsys

bsi.



Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016 & EN ISO 13485:2016

This is to certify that:

JN Medsys Pte Ltd
217 Henderson Road
#02-08
Singapore 159555

Holds Certificate No:

MD 698251

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 & EN ISO 13485:2016 for the following scope:

The design, development, manufacturing, installation, servicing and distribution of analytical instruments, assay kits and reagents for In Vitro Diagnostic applications

For and on behalf of BSI:

Stewart Brain, Head of Compliance & Risk - Medical Devices

Original Registration Date: 2018-11-02
Latest Revision Date: 2018-11-02

Effective Date: 2018-11-02
Expiry Date: 2021-11-01

Page: 1 of 1



...making excellence a habit."

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract. An electronic certificate can be authenticated online. Printed copies can be validated at www.bsi.com/ClientDirectory or telephone +44(0)207007777. Further clarifications regarding the scope of this certificate and the applicability of ISO 13485:2016 & EN ISO 13485:2016 requirements may be obtained by consulting the organization. This certificate is valid only if provided original copies are in complete set.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: +44 345 080 9000
BSI Assurance UK Limited, registered in England under number 2803321 at 389 Chiswick High Road, London W4 4AL, UK.
A Member of the BSI Group of Companies.

PROTECT™ COVID-19 PCR KIT

The **ProTect™ COVID-19 RT-qPCR kit** from JN Medsys is a nucleic acid test that detects for the virus that causes COVID-19. It identifies the genetic material of the SARS-CoV-2 virus isolated from infected patient samples using the reverse transcription quantitative polymerase chain reaction method. Known as the gold standard method, it provides high sensitivity and specificity to identify the virus even at low copies. It allows the accurate identification of the virus isolated from respiratory specimens of infected patients.




Target multiple genes on the virus

The kit is developed based on the United States Center for Disease Control and Prevention (CDC)'s recommended protocol, and detects for multiple genes on the virus RNA, which includes the specific regions of the COVID-19 virus (N1 and N2). Different countries have different protocols for detecting the virus. The reason for targeting multiple genes is to ensure that first, the kit can cover all the different strains of the virus. Second, when the virus starts to mutate over time, and the mutation site unknown, it may be more accurate to detect for multiple targets. The detection of multiple genes may help to increase the accuracy of the tests by reducing the incidences of false negative.

Country	Institute	Gene targets	Number of targets
US	US CDC	Three targets in N gene	3
Germany	Charité	RdRP, E, N	3
China	China CDC	ORF1ab and N	2
France	Institut Pasteur, Paris	Two targets in RdRP	2
Hong Kong SAR	HKU	ORF1b-nsp14, N	2
Japan	National Institute of Infectious Diseases	Pancorona and multiple targets, Spike protein	2
Thailand	National Institute of Health	N	1

Source: WHO

Detection workflow




1. Specimen collection

A nasopharyngeal specimen is collected from the patient using a nasal swab




2. RNA extraction

The specimen is then processed to extract out the RNA from the SARS-CoV-2 virus. This requires the use of a specialised RNA extraction kit and a centrifuge.



3. Polymerase chain reaction (PCR)

The extracted RNA is then put through a process called polymerase chain reaction. Here specific genetic signatures of the SARS-CoV-2 virus are being identified. If it is present, the process makes many copies of the gene fragments. This process requires a specialised PCR kit and a real-time PCR instrument.



4. Results

If the virus is present, the real-time PCR instrument will detect fluorescent signals coming from the PCR reaction. This will confirm the presence of the virus.

Features of the kit



Submitted for US FDA
Emergency Use Authorisation



Sensitivity down to 10
copies/reaction



Approved by Singapore
HSA



Clinically validated at
Singapore hospital



Detects multiple genes
for better accuracy



Competitively priced to
encourage mass usage



Detects all strains of the virus






Results in under 1.5 hours



Manufactured in Singapore
under ISO13485 standards

Performance



Test principle	One-step RT-PCR (Taqman [®] -based)
Limit of detection (LOD)	10 copies/reaction
Precision	< 2%
Dynamic range	6 log
Specificity	Specific to the SARS-CoV-2 virus based on <i>in silico</i> sequence analysis
Efficiency	> 95%


Regulatory approval

Singapore Health Sciences Authority

European Union CE mark

Philippines FDA

Health Sciences Authority
11 Outram Road, Singapore 169078
Tel: 65 6213 0838 Fax: 65 6213 0749
Website: www.hsa.gov.sg



HSA 600/36/01

25 March 2020

JN Medsys Pte Ltd
217 Henderson Road
#02-08
Singapore 159555

Dear Low Huiyu


RE: STATUS OF SUPPLY OF MEDICAL DEVICES IN SINGAPORE

This letter serves to confirm that the following medical device product(s) have been issued Provisional Authorisation for supply in Singapore and may be exported from Singapore.

No.	Device Proprietary Name	Intended Use
1	ProTect™ COVID-19 RT-qPCR Kit (10024) (10027)	The ProTect™ Coronavirus Disease 2019 (COVID-19) RT-qPCR kit by JN Medsys provides all necessary reagents for the in vitro qualitative detection of SARS-CoV-2 from upper respiratory nasopharyngeal specimens. The test is targeting SARS-CoV-2 N1, N2 and N3 genes and the human RNase P control gene.

Product Owner: JN Medsys Pte Ltd
217 Henderson Road
#02-08
Singapore 159555

Manufacturing Site: JN Medsys Pte Ltd
217 Henderson Road
#02-08
Singapore 159555



Page 1 of 2

Health Products Regulation Group • Blood Services Group • Applied Sciences Group
© The Queen's Hospital of Health, The Singapore Health Services & Services & Services

 **Obelis**
EUROPEAN AUTHORIZED REPRESENTATIVE CENTER



**CERTIFICATE
OF
IVD NOTIFICATION**

Ref. No.: BS 9087-2020
Order No.: BS 9053-2020
Date: 24/04/2020

THIS IS TO CERTIFY THAT, ACCORDING TO THE COUNCIL DIRECTIVE 98/79/EC, OBELIS S.A. (O.E.A.R.C.) PERFORMED ALL NOTIFICATION DUTIES AND RESPONSIBILITIES AS THE EUROPEAN AUTHORIZED REPRESENTATIVE (EC REP) OF:

NAME: JN MEDSYS PTE LTD
ADDRESS: 217 HENDERSON ROAD, #02-08, SINGAPORE 159555

AS STIPULATED AND DEMANDED BY THE FOREMENTIONED DIRECTIVE.

The Manufacturer declares that the IVD devices comply with the Directive including all essential requirements.

The Manufacturer has provided Obelis s.a. (O.E.A.R.C.) with all the appropriate declarations according to the 98/79/EC Directive – article 10 requirements including the EC Declaration of Conformity confirming that his In-Vitro Diagnostics medical devices, as stipulated here above, are fulfilling the applicable requirements of the European Council Directive 98/79/EC.

The notification of the following In-Vitro Diagnostic medical devices has been completed by Obelis s.a. (O.E.A.R.C.) on the 23/04/2020 in compliance with the European Council Directive 98/79/EC - article 10 requirements.

IN-VITRO DIAGNOSTIC MEDICAL DEVICES: PLEASE SEE ANNEX A - LIST OF DEVICES (1 PAGE, 1 DEVICE)

As of the 24/04/2020, and as long as the manufacturer will continue complying with the hereabove mentioned requirements* he therefore:

- Is required to affix the CE marking on these devices;
- Place these devices in the Territory of Belgium and/or the other EEA Member States (excluding territories not in alignment with Decision 2010/227/EU).


Mr. G. Elkayam CEO
Obelis sa


   Obelis European Authorized Representative Center is a member of the European Association of Authorized Representatives (E.A.A.R.), ISO 9001 : 2015 and ISO 13485 : 2016 certified in accordance to the profession of a European Authorized Representative.

* This Certificate will be automatically void if the notification is rejected by the EU Authorities or upon termination of the EAA agreement.

Registered Address: BS, General Wafes 53-1000 Brussels / Registered Office Address: B4 Grand Whitlock 30, B-1200 Brussels - Belgium
T: +32 (0) 2 732 5064 | F: +32 (0) 2 732 6009 | Email: mail@obelis.net | Website: www.obelis.net
V2 - C: 0064716 - 32/05/2019





 Republic of the Philippines
Department of Health
FOOD AND DRUG ADMINISTRATION



NO.	PRODUCT NAME	MANUFACTURER
12	2019-nCoV NUCLEIC ACID-BASED DIAGNOSTIC REAGENT KIT (FLOURESCENT PCR)	Sansure Biotech Inc. - No. 680, Lusong Road, Hi-Tech Development Zone, Yuelu District, Changsha, Hunan Province, 410205, People's Republic of China
13	Xpert® Xpress SARS-CoV-2	Cepheid - 904 Caribbean Drive, Sunnyvale, CA 94089-1189 USA
14	NUCLEIC ACID REAGENT TEST KITS FOR NOVEL CORONAVIRUS 2019-NCOV (FLUOROMETRIC PCR)	Shanghai ZJ Bio-Tech Co., Ltd. - Floor 1, Building B, building 20, Floor 1, Building A, building 21, No. 528, Ruigang Road, Zhangjiang high-tech industry east district, Shanghai, China
15	abTES™ COVID-19 REAL-TIME qPCR I KIT	AITbiotech Pte Ltd - 25 Pandan Crescent #05-15 TIC Tech Centre, Singapore
16	SARS-CoV-2 Fluorescent PCR	Maccura Biotechnology Co., Ltd. - 16#, Baichuan Road, High-Tech Zone, 611731 Chengdu, People's Republic of China
17	REAL-TIME FLUORESCENT RT-PCR KIT FOR DETECTING 2019-nCoV	BGI Biotechnology (Wuhan) Co., Ltd. - Research and Development Building B2, Wuhan National Bio-Industry Base B, C, D Suite, No. 666, High-Technology Road, Wuhan Donghu High-Technology Development Zone, People's Republic of China
18	TaqPath COVID-19 COMBO KIT	Thermo Fisher Scientific, Inc. - 5781 Ban Allen Way, Carlsbad, CA 92008, USA
19	ProTect™ COVID-19 RT-qPCR Kit (10024) (10027)	JN Medsys Pte Ltd., -217 Henderson Road #02-08, Singapore
20	cobas® SARS-CoV-2 MAT CODE: 09175431190 TO INCLUDE: cobas® SARS-CoV-2 Control Kit P/N: 09175440190 AND cobas® 6800/8800 Buffer Negative Control Kit P/N: 07002238190	Roche Molecular Systems, Inc. - 4300 Hacienda Drive Pleasanton, CA, USA

* Updated list in yellow highlights

Civic Drive, Filinvest City, Alabang 1781 Muntinlupa, Philippines
Trunk Line +63 2 551 1900 Fax: +63 2 507 0751
Website: www.fda.gov.ph Email: info@fda.gov.ph

Clinical samples validation

COMPARISON STUDIES AT A SINGAPORE HOSPITAL

The ProTect™ kit was validated with clinical samples from a local hospital in Singapore. This hospital has been doing COVID-19 testing using another PCR kit. Overall there is a high degree of agreement (~ 95%) between the two different kits.

		ProTect™ Kit		
		Positive	Negative	Total
Clinical Data From Local Hospital	Positive	27	2	29
	Negative	1	28	29
	Total	28	30	
		Positive Percentage Agreement = 96.4%	Negative Percentage Agreement = 93.3%	
		Overall Rate of Agreement = 94.8%		

ProTect in the news


ASIANSCIENTIST

TOP NEWS | IN THE LAB | HEALTH | TECHNOLOGY | PHARMA | ACADEMIA | FEATURES

Developing Diagnostics In Double-Quick Time

JN Medsys has developed two new coronavirus diagnostic kits to meet the surge in demand for COVID-19 diagnostics.

SHARE
f SHARE
t TWEET
s SHARE



AsianScientist (Mar. 20, 2020) – All over the world, countries are sinking into a public health crisis as the COVID-19 outbreak develops into a global pandemic. Highly infectious, the coronavirus SARS-CoV-2 can spread from person to person even before symptoms appear. Accurate, widespread diagnosis is therefore critical to curbing new viral outbreaks, allowing infected individuals to be identified and quickly quarantined to stop the train of transmission.

South Korea, one of the first countries to experience the epidemic outside China, has shown that widespread testing might be a good solution to contain the virus. The country's elaborate testing program, which saw more than 270,000 people tested thus far, has helped to slow down the influx of new cases since its record high of 909 cases on February 29, 2020.

THE BUSINESS TIMES

551.00 | A SINGAPORE PRESS HOLDINGS PUBLICATION | www.business-times.com.sg | @btnews | CO REGD NO 198402384E | MCI (P) 039/08/2015 | Tuesday April 14, 2020

Camtech, JN Medsys to boost supply of Covid-19 test kits

By Vivienne Tay
vtay@sph.com.sg
@VivienneTayBT

Singapore

TWO Singapore firms, Camtech Diagnostics and JN Medsys, are planning to boost the production of two Covid-19 test kits after obtaining provisional authorisation from Singapore's Health Sciences Authority (HSA).

This comes as demand for diagnostic test kits rises amid the novel coronavirus outbreak in various countries globally.

To meet the projected demand, Camtech Holdings is planning to expand manufacturing in Singapore for the two test kits, according to Kuok Meng-Han, the group's founder and managing director.

Camtech Diagnostics is the research and product development arm of local technology company Camtech Holdings. The latter is an investor in life science firm JN Medsys.

The companies jointly announced on Tuesday that they have obtained provisional authorisation from Singapore's HSA for two Covid-19 test kits – a rapid test kit which produces results in 10 minutes, and a nucleic acid test kit which takes two hours to produce results but has more than 95 per cent accuracy in identifying the coronavirus.

Having a provisional authorisation from HSA allows companies to supply their Covid-19 test kits to healthcare institutions, private hospitals, medical clinics and clinical laboratories in Singapore.

In addition to HSA's provisional licence, JN Medsys has also obtained approval for its test kit from the Food and Drug Administration of the Philippines. It is also in the process of getting approval from the US Food and Drug Administration.

The announcement added that JN Medsys kits have been clinically validated in Singapore and abroad.

Camtech Diagnostics' Covid-19 test is called the Covid-19 IgG/IgM rapid qualitative screening test, which can be a tool in helping to screen and identify asymptomatic carriers of the coronavirus. That said, the test should not be used for confirmatory testing or as a sole basis of diagnosis, the statement noted.

Meanwhile, JN Medsys' test kit is a nucleic acid test that detects the virus which causes Covid-19. Called the ProTect Covid-19 RT-qPCR Kit, the PCR kit identifies the genetic material of the Covid-19 virus isolated from infected patient samples using the reverse transcription quantitative polymerase chain reaction method.

This method is also known as the gold standard method due to its high sensitivity and specificity to identify the virus even at low copies.

The PCR kit also uses protocol from the United States Centers for Disease Control and Prevention.

Johnson Ng, JN Medsys chief executive officer, added that PCR kits have been deployed to China, Malaysia and Indonesia. The PCR kits will be distributed to Asean, Europe and South America soon.

"Digital PCR is already a core product for our existing business, and we have now extended our business line to PCR kits to cater to worldwide demand," he added.

Both Camtech Diagnostics and JN Medsys previously received Spring Singapore innovation grant funding, as well as private investment from local angel investors.

FAQs

- **The kit detects for COVID-19. How does it work?**

The ProTest™ kit detects for the presence of the SARS-CoV-2 virus which causes COVID-19. It does that by detecting for the genetic material of the virus, which in this case is the ribonucleic acid (RNA). The kit targets multiple genes of the virus. Once it “finds” these genes, it makes many copies of that genes using a process known as polymerase chain reaction. As more copies of those genes are being made, it generates fluorescent signals that can be detected using a real-time PCR instrument. The PCR is a very accurate and sensitive test as it can detect for very low copies of the virus. Currently, the detection of the SARS-CoV-2 virus using PCR is the recommended method used by most countries for confirmation of COVID-19 infection in patients.

- **What is the entire process for detecting for the virus?**

First, a nasopharyngeal specimen is collected using a nasal swab. The swab is then immersed into a transport medium. Second, RNA extraction is carried out on the specimen. This process lyses the viral particles to release the RNA and then purifies it. This step must be carried out with an RNA extraction kit and a bench-top centrifuge. Third, the extracted RNA is put through the PCR process. If specimen contains the SARS-CoV-2 virus, the real-time PCR process will show a positive result, confirming the presence of the virus.

- **How long does it take?**

The entire process from RNA extraction to PCR can be conducted in under 2 hours, with the PCR process requiring about 1h 15min.

- **How is this different from the rapid test kits that can be performed in 15 minutes?**

The rapid test kits that can generate results in 15 min detects for the antibodies present in an infected person. An infected person takes about a week to generate these antibodies in response to the viral infection. During this period, the rapid test may return a negative result due to the lack of antibodies even though the patient might be infected.

FAQs

- **What equipment do I need?**

A real-time PCR instrument is needed. We recommend the following:

- ThermoFisher Scientific QuantStudio series (3, 5, 6 and 7) instruments and 7500 real-time system
- Bio-Rad CFX96 Touch™ Real-time PCR System
- Roche LightCycler® 96 and 480 Systems

- **On which real-time PCR instrument has this been validated? Is it compatible with others?**

The kit has been validated on ThermoFisher Scientific's QuantStudio 3 and Roche's LightCycler 480 II.

- **How is this kit different from the Fortitude kit from A*Star?**

Both the ProTect™ and A*Star's Fortitude 2 use the PCR method to detect for the SARS-CoV-2 virus. The Fortitude 2 detects for a single gene on the virus. The ProTect™ detects for multiple genes on the virus, in accordance to the US Communicable Disease Center (CDC) protocol. Detecting for multiple genes on the virus ensures that the kit can detect 100% of all the different strains of the virus. It also ensures that the kit can still detect the virus even if it mutates.

- **Which targets on the SARS-CoV-2 virus does it detect?**

It detects for the N gene of the SAR-CoV-2 virus, in accordance to US CDC recommended protocol.

- **Each kit has how many tests?**

Each kit has 100 tests.

FAQs

- **What is the dimension and weight of each kit?**

Each is contained in 2 small boxes. Each box is 73 x 49 x 54 (H) mm. Two boxes bundled together will be 73 x 98 x 54 (H) mm, weighing 100g

- **What is the storage condition for the kits?**

The kits are to be stored at -20°C and repeated freeze thawing of the reagents are to be avoided.

- **What is the shipping condition for the kits?**

The kits are shipped in dry ice in a temperature controlled carton box. Each small carton measures 54 x 42 x 40 cm and weighs 25kg (including the dry ice) and can fit about 12 kits. For larger volume, we ship in specialised containers in a pallet.

- **How much is it to ship 1,000 tests?**

The price varies depending on destination. Price can also fluctuate depending on availability of flights.

- **What countries can we export to?**

We currently can export to most places worldwide, subject to availability of freight.

- **What should I do if I need a bulk order of a million test?**

Due to the size of the order, please contact us to discuss the delivery schedule

- **What is the expected lead time?**

Lead-time is dependent on the quantity of orders.



The **ProTect™ COVID-19 PCR kit** provides all the necessary reagents for the accurate and sensitive detection of the Sars-CoV-2 virus using reverse transcription real-time polymerase chain reaction (RT-qPCR).

Contact sales@jnmedsys.com to order now or visit <https://jnmedsys.com/kits/> to find out more.



Accurate



Made in Singapore

IVD

**Provisional
authorisation from
Singapore HSA**



Affordable



**Based on US CDC
protocol**